MAGUS CHD50 DIGITAL CAMERA

MAGUS

Levenhuk Inc. (USA):
928 E 124th Ave. Ste D, Tampa, FL 33612, USA,
+1-813-468-3001,
contact_us@levenhuk.com

Levenhuk Optics s.r.o. (Europe):
V Chotejně 700/7, 102 00 Prague 102, Czech Republic,
+420 737-004-919,
sales-info@levenhuk.cz

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EN MAGUS CHD50 DIGITAL CAMERA

Never look at the sun, a bright light source, or laser beam through the camera. It can cause permanent damage to your eyes.

MAGUS high-speed HDMI cameras are designed to be paired with a digital microscope. Sony Starvis high-sensitivity sensor ensures perfect color reproduction and image quality even in low light conditions. You can install the camera software to view, save, and edit captured images and videos. The rugged aluminum housing is designed for operation at low temperatures (−10°C).

HDMI

1. Use an HDMI cable to connect the camera to an HDMI monitor (4).
2. Connect a USB mouse to the connector (7).
3. Insert the memory card into the memory card slot (2). Or insert the USB memory stick into the slot (5).
4. Connect the camera to the power supply (3). Turn the camera on (1).
5. Turn the monitor on and watch videos in the built-in program.

CONNECT THE CAMERA TO A PC/LAPTOP VIA USB2.0

1. Install the MAGUSView software from the USB flash drive (included in the kit) on the PC.
2. Connect the camera to the power supply (3). Turn the camera on (1).
3. Connect the camera to a PC with a cable by inserting one plug into the connector (3) and the other plug into a USB port on the PC.
4. Open the MAGUSView program. Choose a camera from the list.

CONNECT THE CAMERA TO A PC/LAPTOP VIA WI-FI

Make sure your PC supports Wi-Fi.

1. Install the MAGUSView software from the USB flash drive (included in the kit) on the PC.
2. Connect the camera to the power supply (3). Turn the camera on (1).
3. Insert the Wi-Fi antenna to the connector (2).
4. After the antenna starts blinking, select the Wi-Fi network that you want to connect to on your PC. The password is 12345678.
5. Open the MAGUSView program. The camera will be automatically recognized by the camera software.

CONNECT MULTIPLE CAMERAS TO A PC/LAPTOP VIA WI-FI

Make sure your PC supports Wi-Fi.

1. Install the MAGUSView software from the USB flash drive (included in the kit) on the PC.
2. Connect the camera to the power supply (3). Turn the camera on (1).
3. Insert the Wi-Fi antenna to the connector (2).
4. After the antenna starts blinking, select the Wi-Fi network that you want to connect to on your PC. The password is 12345678.
5. To connect multiple cameras:
Move the mouse cursor to the bottom of the window, the control toolbar will automatically pop up. Select Settings>Network> WLAN, then select the WLAN STA mode. Enter the SSID and the password of the connected router.
6. Open the MAGUSView program. The camera will be automatically recognized by the operating system. You can view the images from the multiple connected cameras by choosing the required camera from the list in the left part of the video display window.

CALIBRATION USING A CALIBRATION SLIDE


You should use a special stage micrometer (calibration slide) to determine the linear dimensions (in millimeters or microns).

The calibration slide is a transparent glass (of the same size as the specimen slide) that has a micrometer scale with a scale division of 0.01mm etched on the surface.

The calibration slide (stage micrometer) is used to calibrate the image analysis software for measurements in actual units. In the calibration mode, you should capture an image of the micrometer scale with every objective magnification and indicate the known distance. That lets you establish a scale of the image in actual units (micrometer, millimeter, etc.).

Calibration:


1. Place the calibration slide on the microscope stage.
2. Select the required objective and set the maximum camera resolution.
3. Get a contrast image of the scale on the monitor screen and capture the image.
4. Select the "Calibrate" function in the software you are using.

- 5. Double-click on the maximum visible distance and enter the value in actual units.
- 6. Enter the calibration setting and check the result. The program will save the calibration factor.
- 7. You can select any measurement unit later, and all the results will be re-calculated in accordance with this selection.


CONTROL PANEL


To open the Camera control panel, move the mouse cursor to the left of the window. The control panel will automatically pop up.

CAMERA CONTROL PANEL	FUNCTION	FUNCTION DESCRIPTION
	Snap	Captures the image and saves it to the SD card.
	Record	Records a video and saves it to the SD card.
	Auto Exposure	Automatically adjusts the exposure time and gain according to the exposure compensation value.
	Exposure Compensation	Available when Auto Exposure is selected. Adjusts the exposure compensation according to the current video brightness to achieve the proper brightness value.
	Exposure Time	Available when Auto Exposure is not selected. Adjusts the time of exposure and thus adjusts the brightness of the video.
	Gain	Adjusts the brightness of the video. The Noise will be adjusted accordingly.
	Red	Adjusts the proportion of the red color in RGB in the video.
	Green	Adjusts the proportion of the green color in RGB in the video.
	Blue	Adjusts the proportion of the blue color in RGB in the video.
	Auto White Balance	Adjusts the White Balance of the video continuously.
	Manual White Balance	Adjusts the Red color or Blue color parameters to set the White Balance of the video. White balance is adjusted to obtain a natural ratio of cold and warm colors on the image with regard to the light source type. Move the corresponding sliders to adjust the values.
	ROI White Balance	Adjusts the white balance of the Region of Interest when the Region of Interest content is changed.
	Sharpness	Adjusts the Sharpness level of the video.
	Denoise	Removes the noise from the video.
	Saturation	Adjusts the Saturation level of the video.
	Gamma	Allows you to adjust the midtone values of the image. The higher the Gamma value is, the darker the image will be.
	Contrast	Adjusts the Contrast level of the video.
	DC	Allows you to use a light source without light flickering.
	AC(50HZ)	Eliminates light flickering when a 50Hz light source is used.
	AC(60HZ)	Eliminates light flickering when a 60Hz light source is used.
	Default	Restores the default settings.

MEASUREMENT TOOLBAR

To open the Measurement Toolbar, move the mouse cursor to the top of the window. The toolbar will automatically pop up.

ICON	FUNCTION	ICON	FUNCTION
	Float/Fix switch of the Measurement Toolbar		Circle: center and radius
	Show/Hide Measurement Objects		3 Point Circle: drawing up of a circle using three points
	Select the desired Measurement Unit		Annulus
	Select Magnification for Measurement after Calibration		3 Points Annulus
	Object Select		Two Circles and its Center Distance: drawing up of two circles and the measurement of the distance between the centers of these two circles
	Angle		Two 3 Point Circles and the Center Distance: drawing up of two three-point-circles and the measurement of the distance between the centers of these circles
	4 Point Angle		Arc
	Point		Text
	Arbitrary Line		Polygon
	3 Point Line		Curve
	Horizontal Line		Scale Bar
	Vertical Line		Arrow
	Parallel		Executes Calibration with the help of a micrometer to determine the corresponding relation between magnification and resolution as well as between the measurement unit and the sensor pixel size.
	3 Point Vertical Line		Exports the Measurement information to CSV file (*.csv)
	Rectangle		Measurement Setup
	3 Point Rectangle		Deletes all of the measurement objects
	Ellipse		Exits from Measurement mode
	5 Point Ellipse: drawing up of an ellipse using five points	      Use these icons to select Move Left, Move Right, Move Up, Move Down, Color Adjustment, or Delete.	

Select a figure to draw up; click at the desired initial point on the image, click at the desired end point. The selected figure will be drawn up on the image and the figure parameters will be displayed.


CAMERA CONTROL TOOLBAR

To open the Camera Control Toolbar, move the mouse cursor to the bottom of the window. The control toolbar will automatically pop up.


ICON	FUNCTION	ICON	FUNCTION
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color/Gray		Video Freeze
	Display Cross Line		Image Overlay
	Auto Focus		Compare the image with the current video
	Browses the images and videos on the SD Card		Settings
	Check the Version of MAGUSView		

AUTO FOCUS CONTROL PANEL

Auto Focus	When the Auto Focus function is selected, the system will start autofocusing on the object.
Manual Focus	When the Manual Focus function is selected, you need to adjust the camera sensor position manually by using the mouse to scroll the mouse wheel up and down until the object is in focus.
One Push AF	Click the One Push AF button to perform the autofocusing operation just once.
Conjugated Calibration	Click the Conjugated Calibration button to reset the camera sensor to the standard C-mount position. The Conjugated Calibration function helps to calibrate the camera sensor position when the camera video display window is clear and the view from microscope eyepiece is clear. Apply the Conjugated Calibration function when starting the camera for the first time to set the camera sensor at the standard C-mount position. It is required to set the object plane, the eyepiece image plane and the camera adapter image plane at the standard positions. If the height of the object changes, you must make sure that the camera sensor is at the standard C-mount position. Adjust the coarse and fine focus of the microscope to focus on the object. Prior to perform measurements, please start the Conjugated Calibration function to ensure the accuracy of the measurement results (see the Conjugated Calibration option in the Measurement Toolbar Section).


The Focus Region is used to select the region of interest (ROI) for the autofocus operation. You can click any part of the video display window to change the focus region position for the Auto Focus operation.

When you close the Auto Focus Control Panel, the Focus Region will also be closed automatically.

If the Focus Region function is enabled, the Measurement Toolbar will not appear when the mouse cursor is moved to top of the video display window.


SETTINGS

NETWORK

Select the AP mode to connect one digital camera using the Wi-Fi adapter (included).

Select the STA mode to connect multiple digital cameras using the Wi-Fi router. Enter the SSID of the router and the password of the router to start operation.

MEASUREMENT

Global	Precision	Sets the required number of digits after the decimal point. The measurement result value will be displayed with the required precision.
Calibration	Line Width	Defines width of the lines for calibration.
	Color	Defines color of the lines for calibration.
	EndPoint	Type: defines shape of the endpoints of lines for calibration: Null means no endpoints, Rectangle means rectangle shape of endpoints. It facilitates the object alignment.

Object Select, Angle, 4 Points Angle, Point, Arbitrary Line, 3 Points Line, Horizontal Line, Vertical Line, Parallel, 3 Points Vertical Line, Rectangle, 3 Point Rectangle, Ellipse, 5 Point Ellipse, Circle, 3 Points Circle, Annulus, 3 Points Annulus, Two Circles and Its Center Distance, 3 Points Two Circles and Its Center Distance, Arc, Text, Polygon, Curve, Scale Bar, Arrow

You can Left-click the  and select Measurement to open the corresponding window of settings. You can specify the parameters of the required measurement objects.

MAGNIFICATION

Name	Sets the magnification item name when calibration, such as 4x, 10x, 100x, etc. Other information could be added too, for example, microscope model, operator name, etc.
Resolution	Sets the pixels per meter.
Clear All	Clears all of the calibrated magnification items.
Delete	Deletes the selected magnification item.
Up	Click Up to move up the selected magnification item up.
Down	Click Down to move the selected magnification item down.

IMAGE FORMAT

Image Format	JPEG: you can save files in the JPEG format to get a high compression and high quality images and to save the storage space. TIFF: Tag Image File Format (TIFF) is a flexible bitmap format that is used to store images including photos and artistic images.
Measurement Object Saving Method	Burn-in Mode: merges the measurement objects into one current image. Further edition of the measurement objects is impossible. The changes are irreversible. Layered Mode: saves the measurement objects at different layers of the current image into the target file. The measurement objects in the target file can be edited with the use of corresponding software on the PC. The changes are reversible.

VIDEO

Playback Fast Forward/ Reverse Interval	Sets the time interval of the playback of video files.
Video Encode	You can select the video encoding format H264 or H265.

STORAGE

Storage Device	You can save data on the SD card or on the USB flash drive.
File System Format of the Storage Device	You can select the file system format: FAT32, exFAT or NTFS. Use PC to format the SD cards and switch between the file system formats. The maximum video file size of a single file is 4GB. Unknown Status: SD card or USB flash drive not detected or the file system is not identified. FAT32: The maximum video file size is 4GB. exFAT: The maximum video file size is 16EB. NTFS: The maximum video file size is 2TB.

For USB flash drive, it is recommended to use the USB 3.0 interface.


FILES

Image File Name	Auto: Saves the image file automatically with the specified prefix. Manual: You need to name the image file before saving.
Video File Name	Auto: Saves the video file automatically with the specified prefix. Manual: You need to name the video file before saving.

LANGUAGE

You can also select a language from the list.

MISCELLANEOUS

Clarity factor	Displays or hides the clarity factor (sharpness factor) in the video display window, to check the camera focusing.
Ruler	Displays or hides the Ruler in the video window.
Measurement	Displays or hides the Measurement toolbar in the video window.
Overlay	Displays or hides the saving graphics Overlay information in the Burn in Mode.
Grids	Displays or hides a mesh in the fusion mode.

USB video output switch back to mouse operation	Select Automatic restart or Manual restart to use the USB video interface for mouse operation.
ROI Color	Chooses the line color for a rectangle ROI (region of interest)
Cursor	Changes the cursor size according to the screen resolution or personal preferences.
Auto Exposure	Defines the maximum exposure time for automatic settings.
Auto Exposure Region	Select the Auto Exposure reference area.
Camera Parameters Import	Imports the camera parameters from the SD card or from the USB flash drive to use the previously exported camera parameters.
Camera Parameters Export	Exports the camera parameters to the SD card or to the USB flash drive to use the previously exported camera parameters.
Reset to factory defaults	Restores the camera parameters to the default settings.

THE KIT INCLUDES

MAGUS digital camera, 12V, 1A power adapter, HDMI cable, USB mouse, 32GB SD card, USB3.0 cable, USB flash drive with drivers and software, USB3.0 cable, CF-812AC USB Wi-Fi adapter, mounting plate and screws, user manual, and a warranty card.

SYSTEM REQUIREMENTS

Windows 8/10/11 (32bit and 64bit), Mac OS X, Linux, up to 2.8GHz Intel Core 2 or higher, minimum 4GB RAM, USB2.0 port, 19" or larger display.

SPECIFICATIONS

	MAGUS CHD50
Maximum resolution	3840x2160px
Megapixels	8
Sensor	1/1.8" (7.68x 4.32mm) Sony CMOS sensor
Pixel size	2.0x2.0µm
Light sensitivity	505mv with 1/30s
Exposure	0.04–1000ms
Video recording	+
Frame rate, fps at resolution	30fps@3840x2160px (HDMI) 30fps@3840x2160px (WLAN) 30fps@3840x2160px (USB)
Image format	*.jpeg, *.tiff
Video format	*.mp4
Shutter type	ERS
Auto focus	+
Interface	HDMI, Wi-Fi, USB2.0
USB Wi-Fi adapter	802.11n 150 Mbps
Body	metal
Operating temperature range	–10...50°C

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

CARE AND MAINTENANCE

- Do not disassemble the camera on your own.
- Keep the camera away from moisture; do not use it in rainy weather.
- Keep the camera away from shock or excessive pressure.
- Do not overtighten the locking screws.
- Keep the camera away from hazardous environment, home and car heaters, incandescent lamps or open fire.
- When cleaning the lenses, first blow any dust or debris off the surface or wipe them off with a soft brush. Then wipe the lens with a soft cleaning cloth slightly damp with spirits or ether.
- Seek medical advice immediately if a small part or a battery is swallowed.

MAGUS WARRANTY

Magus optic products carry a **5-year warranty** against defects in materials and workmanship. All Magus accessories are warranted to be free of defects in materials and workmanship for **six months** from the purchase date. The warranty entitles you to the free repair or replacement of the Magus product in any country where a Levenhuk office is located if all the warranty conditions are met.

For further details, please visit: www.magusmicro.com

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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MAGUS® is a registered trademark of Levenhuk, Inc.
© 2006–2024 Levenhuk, Inc. All rights reserved.