Levenhuk biological microscopes

Levenhuk 625 (binocular)

Levenhuk 670T (trinocular)

Levenhuk D670T (trinocular)

User Manual

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

levenhuk[®]
Zoom&Joy

- 1 Eyepieces
- 2 Binocular/Trinocular head
- 3 Revolving nosepiece
- 4 Locking screw
- 5 Objective lenses
- 6 Arm
- 7 Stage
- 8 Specimen holder
- 9 Condenser adjustment knob
- 10 Focus adjustment knob
- 11 Abbe condenser
- 12 Coaxial stage control
- 13 Base

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General Information

Levenhuk 625, Levenhuk 670T and Levenhuk D670T biological microscopes are modern state of the art optical instruments made with great attention to detail. The binocular head of 625 and trinocular heads of 670T and D670T are inclined at a 30° angle to the stage that moves up and down with precise coarse and fine focus systems.

These instruments are shipped with high-quality achromatic objective lenses, wide FOV eyepieces and a variable brightness illumination source that together provide for the best quality of resulting images.

With their great quality and reasonable price Levenhuk 625, Levenhuk 670T and Levenhuk D670T are a perfect choice for specialists from all scientific fields. They perform well during clinical research and demonstrations, in universities and laboratories, and may be used for scientific research in agriculture and microbiology.

Specifications

	Levenhuk 625	Levenhuk 670T	Levenhuk D670T
Type	Biological	Biological	Digital
Head	Binocular, inclined at 30°	Trinocular, inclined at 30°	Trinocular, inclined at 30°
Revolving nosepiece	Yes	Yes	Yes
Magnification	40x-2000x	40x-2000x	40x-2000x
Objective lenses	achromatic 4x, 10x, 40x (S), 100x (S, oil immersion)	achromatic 4x, 10x, 40x (S), 100x (S, oil immersion)	achromatic 4x, 10x, 40x (S), 100x (S, oil immersion)
Eyepieces	WF10x / ø 18 mm or 20 mm, WF20x / ø ≥ 11 mm	WF10x / ø 18 mm or 20 mm, WF20x / ø ≥ 11 mm	WF10x / ø 18 mm or 20 mm, WF20x / ø ≥ 11 mm
Tubus length	160 mm (6.3 in)	160 mm (6.3 in)	160 mm (6.3 in)
Stage size	125 x 115 mm (4.92 x 4.53 in)	125 mm x 115 mm (4.92 x 4.53 in)	125 mm x 115 mm (4.92 x 4.53 in)
Stage movement range	60 mm x 30 mm (2.36 x 1.18 in)	60 mm x 30 mm (2.36 x 1.18 in)	60 mm x 30 mm (2.36 x 1.18 in)
Condenser	Abbe with 1.25 N.A. iris diaphragm and a filter	Abbe with 1.25 N.A. iris diaphragm and a filter	Abbe with 1.25 N.A. iris diaphragm and a filter
Focus system	coaxial, coarse and fine	coaxial, coarse and fine	coaxial, coarse and fine
Illumination	6V/20W halogen lamp (variable brightness)	6V/20W halogen lamp (variable brightness)	6V/20W halogen lamp (variable brightness)
Weight	8 kg (17.64 lbs)	4 kg (8.82 lbs)	4 kg (8.82 lbs)
Dimensions (HxWxD)	360 x 160 x 260 mm (14.17 x 6.3 x 10.24 in)	375 x 174 x 220 mm (14.76 x 6.85 x 8.66 in)	375 x 174 x 220 mm (14.76 x 6.85 x 8.66 in)
Power source	separate cable	separate cable	separate cable

The kit

Levenhuk 625: microscope body, binocular head, achromatic objective lenses (4x, 10x, 40x S, 100x S (oil immersion)), WF10x and WF20x eyepieces (2 pcs), 6V/20W halogen lamp, blue filter, vial of immersion oil, dust cover.

Levenhuk 670T: Microscope body, trinocular head, achromatic objective lenses (4x, 10x, 40x S, 100x S (oil immersion)), WF10x and WF20x eyepieces (2 pcs), 6V/20W halogen lamp, blue filter, vial of immersion oil, dust cover.

Levenhuk D670T: Microscope body, trinocular head, achromatic objective lenses (4x, 10x, 40x S, 100x S (oil immersion)), WF10x and WF20x eyepieces (2 pcs), 6V/20W halogen lamp, blue filter, vial of immersion oil, 5 mpx digital camera, software & drivers CD, dust cover.

Parts of the microscope


Levenhuk 625 / Levenhuk 670T / Levenhuk D670T is comprised of nine parts:

- **Base:** not only does it hold the weight of the microscope, it also houses the illumination source, electronics and control mechanisms.
- **Arm:** this piece holds the base, the stage and the head of the microscope together. Coarse and fine focus systems provide for smooth vertical movements of the stage.
- **Rack-and-pinion mechanism:** mounted on the arm, the stage with the condenser are moving vertically along this column. For additional precision, a condenser may be adjusted separately.
- **Head:** a binocular (Levenhuk 625) or a trinocular (Levenhuk 670T / Levenhuk D670T) head is mounted at a 30° angle at the upper end of the arm.
- **Eyepieces:** wide FOV eyepieces WF10x and WF20x are used in these microscopes.
- **Revolving nosepiece:** quadruple revolving nosepiece allows you to change objective lenses smoothly and easily.
- **Objective lenses:** high-quality achromatic objective lenses with 4x, 10x, 40x and 100x magnifications provide for sharp and bright images.
- **Stage:** sturdy and reliable stage with a specimen holder that can be used to move your slides while observing them.
- **Condenser:** Abbe condenser, with 1.25 N.A. iris diaphragm.

Operating principle and illumination

1. **Image creation system:** objective lens (6), prism (7) and eyepiece (9). The objective lens (6) magnifies a specimen (0), light rays pass through a prism (7), refract at a 45° angle and create an image in the eyepiece. Total magnification may be calculated by multiplying magnifications of the eyepiece and the objective lens used.

2. **Illumination system:** lamp (1), collector lens (2), diaphragm (3) and condenser (4). Light emitted from a lamp (1) passes through a collector lens (2) and illuminates a diaphragm (3). After this, it is focused by a condenser (4). This illumination system is used for observations of a specimen (0) in transmitted light. However, you can also use a different type of illumination for observations in reflected light.

Digital camera

Levenhuk D670T comes with a 5 mpx C510T NG digital camera.

The camera allows you to observe specimens in fine detail and true colors on your PC monitor and save images on the hard drive.

The special software that comes in the kit allows you to view and edit the resulting images.

Supported file formats include: *.bmp, *.jpg, *.jpeg, *.png, *.tif, *.tiff, *.gif, *.psd, *.ico, *.emf, *.wmf, etc.

Connect the camera to your PC via a USB-cable that also acts as a power cable, so you do not have to plug it into a power source.

Compatible OS: Windows 2000/XP/2003/Vista

Camera model	C510 NG
Max. resolution (for still images)	2592x1944
Megapixels	5.0
Sensor	1/2.2" CMOS
Dimensions (cylinder unit maximum diameter), mm	50
Pixel size	2.2µm x 2.2µm
Sensitivity, V/lux.sec@550 nm	1.2
Mounting location	Eyepiece tube (replaces the eyepiece)
Image format	BMP, TIFF, JPG, PICT, SFTL, etc.
Field of view diameter	18 mm
Spectral range	400 nm - 650 nm
Exposure	ERS
White balance	Auto/manual
Exposure control	Auto/manual
Software image control	Image size, brightness, exposure control
Interface	USB 2.0, 480 Mb/s
System requirements	Windows 2000/XP(SP2)/2003/Vista; USB port
Software	USB 2.0 driver, software for image viewing and editing
Shell	Metal cylinder, painted black
Power source	USB 2.0 cable
Temperature range	-30... 70 °C (-22... 158 F)

Caution

- Never look at the sources of bright light or lasers through the camera or the microscope: it will cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS!
- Do not disassemble the camera on your own.
Protect the camera from moisture; do not use it in rainy weather.
- Protect the camera from impacts or excessive pressure.
Do not overtighten the locking screws.
- Keep the camera away from hazardous environment, home and car heaters, incandescent lamps or open fire.
- When cleaning the lenses, first blow any dust or debris off the surface or wipe them off with a soft brush. Then wipe the lens with a soft cleaning cloth dipped in lens cleaning fluid.
- Never touch the optics with your fingers.

Care and maintenance

After unpacking your microscope and before using it for the first time check for durability of every component and connection. Do not apply excessive pressure, lest you break your instrument. Be careful during your observations; always replace the dust cover after you are finished with observations to protect your microscope from dust and stains. Do not try to take your microscope apart, as you might damage the intricate optical system within.

Store the microscope in a dry and cool place, away from possible impacts and dust. Do not try to perform repairs of any kind on this instrument - contact a Levenhuk service center instead. If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope itself.

Using the microscope

- Insert eyepieces into the binocular head and attach objective lenses to the revolving nosepiece in ascending order. Place a specimen in the middle of the stage and hold it in place with stage clips.
 - Connect your microscope to a power outlet and adjust the brightness of illumination. When you finish your observations, always dim the illumination before unplugging the instrument. Start at the lowest magnification and center the observed specimen. Afterwards, switch to a
 - higher magnification by rotating the revolving nosepiece. If the image is out of focus, you can use the fine focus system to adjust it. When using a 100x objective lens, you have to fill the space between the lens and the specimen with immersion oil. There should be no air bubbles in
 - the oil, as this will reduce the resulting quality of the image.
- For the best quality of your views, illumination has to be carefully adjusted. When switching to a different objective lens, you will most likely have to readjust the iris diaphragm and the illumination brightness.
- If you need to change the light bulb, unplug your microscope, wait for the bulb to cool and change it. Note that the lamp has to be leveled and inserted all the way to the end.

International warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a lifetime warranty against defects in materials and workmanship. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for six months from date of retail purchase. Levenhuk will repair or replace such product or part thereof which, upon inspection by Levenhuk, is found to be defective in materials or workmanship. As a condition to the obligation of Levenhuk to repair or replace such product, the product must be returned to Levenhuk together with proof of purchase satisfactory to Levenhuk.

This warranty does not cover consumable parts, such as batteries.

A Return Authorization (RA) Number must be obtained in advance of return. Contact the local Levenhuk branch to receive the RA number to be displayed on the outside of your shipping container. All returns must be accompanied by a written statement setting forth the name, address and telephone number of the owner, including a description of any claimed defects. Parts or products for which replacement is made will become the property of Levenhuk.

The customer will be responsible for all costs of transportation and insurance to and from Levenhuk or its authorized dealers and will be required to prepay such costs.

Levenhuk will use reasonable efforts to repair or replace any product covered by this warranty within thirty days of receipt. If a repair or replacement will require more than thirty days, Levenhuk will notify the customer accordingly. Levenhuk reserves the right to replace any product that has been discontinued from its product line with a new product of comparable value and function.

This warranty does not apply to any defects or damages resulting from alteration, modification, neglect, misuse, usage of improper power sources, damage in transportation, abuse, or any cause other than normal use, or to malfunction or deterioration due to normal wear. Levenhuk disclaims all warranties, express or implied, whether of merchantability or fitness for a particular use, except as expressly set forth herein. The sole obligation of Levenhuk under this limited warranty will be to repair or replace the covered product, in accordance with the terms set forth herein. Levenhuk disclaims liability for any loss of profits, loss of information, or for any general, special, direct, indirect or consequential damages which may result from breach of any warranty, or arising out of the use or inability to use any Levenhuk product. Any warranties which are implied and which cannot be disclaimed will be limited in duration to a term of six months for accessories from the date of retail purchase.

Some states/provinces do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations and exclusions may not apply to you. This warranty gives you specific legal rights, and you may have other rights which vary from state to state or province to province.

Levenhuk reserves the right to modify or discontinue any product without prior notice.

NOTE: This warranty is valid to USA and Canadian customers who have purchased this product from an authorized Levenhuk dealer in the USA or Canada. Warranty outside the USA or Canada is valid only to customers who purchased from an authorized Levenhuk dealer in the specific country or international distributor. Please contact them for any warranty service.

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch:

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Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2006-2014 Levenhuk, Inc. All rights reserved.
Levenhuk, Inc. 1935 Brandon Court, Suite A-1