TR

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Microscope Digital Camera M PLUS SERIES

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
M500 PLUS
M800 PLUS
M1000 PLUS
M1200 PLUS
M2000 PLUS

www.levenhuk.com

Levenhuk Inc. (USA): 928 E 124th Ave. Ste D, Tampa, FL 33612, USA, +1 813 468-3001, contact_us@levenhuk.com
Levenhuk Optics s.r.o. (Europe): V Chotejně 700/7, 102 00 Prague 102, Czech Republic, +420 737-004-919, sales-info@levenhuk.cz
Levenhuk® is registered trademark of Levenhuk, Inc.
© 2006–2025 Levenhuk, Inc. All rights reserved.
20250529

EN Microscope Eyepiece Digital Cameras Levenhuk M PLUS Series

Never look at the sources of bright light or lasers through the camera or the microscope: it will cause PERMANENT DAMAGE TO YOUR EYES!

Levenhuk eyepiece digital cameras have been designed exclusively for microscopes. They can be installed on any type of optical microscopes, be it biological or metallurgical, mono or stereo microscope. The observed view may be transferred directly to your PC, and with the software included in the kit you can view and edit the resulting image.

Features

- Levenhuk M PLUS Series includes cameras with resolutions ranging from 5 to 20 megapixels;
- The camera allows you to watch images on the screen in a real-time mode, and to capture photos and videos;
- The camera renders realistic colors with high level of detail;
- Levenhuk M PLUS Series cameras are compatible with microscopes of 23.2mm eyepiece tube diameter (the camera is installed in place of the eyepiece);
- The standard kit includes Levenhuk software necessary to view and edit your images (brightness, image size, etc.);
- The camera can be charged and connected to your PC with USB-cable.

Package

Levenhuk digital camera, adapter 23.2 mm, 2 adapters for microscopes with different eyepiece tube diameters (23.2–30mm, 23.2–30.5mm), USB cable, installation CD (setup wizard, drivers, software), user manual and warranty.

System requirements

Mac OS 10.12, Linux Ubuntu 14.04, Windows XP/Vista/7/8/10/11 (32-bit or 64-bit), CPU IntelCore 2 or above; RAM 2GB, USB port 2.0, CD-ROM.

Care and maintenance

- Never look at the sources of bright light or lasers through the camera or the microscope: it will cause PERMANENT DAMAGE TO YOUR EYES!
- Do not disassemble the camera on your own.
- Keep the camera away from moisture; do not use it in rainy weather.
- Keep the camera away from shock or excessive pressure.
- Do not overtighten the locking screws.
- Keep the camera away from hazardous environment, home and car heaters, incandescent lamps or open fire.
- When cleaning the lenses, first blow any dust or debris off the surface or wipe them off with a soft brush. Then wipe the lens with a soft cleaning cloth slightly damp with spirits or ether.
- Seek medical advice immediately if a small part or a battery is swallowed.

Warranty: lifetime. For further details, please visit our web site:

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